



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Disposable Laboratory Coats

**Product Code :**  
LBNY-00017-1100013



**Description :**

Disposable Lab Coats come complete with a knitted collar and cuffs. Designed to protect healthcare professionals, lab technicians and other users from fluids and cross-contamination. Made from high quality, ultra-soft fabric designed to keep you protected and comfortable. Premium protective apparel is available in both lab coats and hip length jackets. High quality front snap button closure  
Comfortable knit collar and cuffs  
Fluid-resistant material meets AAMI L2 standards  
Anti-static coating  
3-pocket design  
Convenient side slits for easy access to pockets.  
42 GSM high quality, breathable SMS material.

**Technical Specification :**

Disposable Laboratory Coats



**Laboratory Instrument India**

