



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Lightweight Cotton Inspection Gloves

**Product Code :**  
LBNY-00017-1100015



**Description :**

Ideal for general inspection work, or for use as liners underneath unlined liquid proof gloves.  
Wear alone for inspection purposes  
100% cotton lisle gloves  
Two-piece construction  
Ambidextrous.  
Can be used as liners under liquid proof gloves.

**Technical Specification :**

Lightweight Cotton Inspection Gloves



**Laboratory Instrument India**