



Email : sales@laboratoryinstrumentindia.com

Product Name :
Silicone Hand Protector

Product Code :
LBNY-00017-1100017



Description :

Ideal for handling extremely hot or cold objects, including dry ice.
Provides protection from temperatures ranging from -50°C to 260°C (-70°F to 500°F).
Silicone rubber conforms to the hand and object being handled.
Gripping surface is studded with small nubs for protection and sure handling of objects with temperatures ranging from -50°C to 260°C.
Gripping surface studded with small nubs for positive grip and temperature protection
Tab for hanging near work station.
One size fits all hands, left or right.

Technical Specification :

Silicone Hand Protector



Laboratory Instrument India