



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Unisex Specialized Cuffed Lab Coat

**Product Code :**  
LBNY-00017-1100021



**Description :**

This lab coat is off the charts with all our best material production features. Featuring a 6-gripper closure, neck gripper, and a chest pocket, this Unisex Specialized Cuffed Lab Coat looks as professional as it is durable. Plus, who doesn't need the storage offered by a chest pocket and two lower pockets? Put this lab coat to the test and see why it's a customer favorite. Lab, medical facility, and clinic workers will appreciate the knit cuffs that keep sleeves from interfering with work and extended back coverage to protect from messy jobs.

**Technical Specification :**

Unisex Specialized Cuffed Lab Coat



**Laboratory Instrument India**