



Email : sales@laboratoryinstrumentindia.com

Product Name :
Professional Nitrile Gloves

Product Code :
LBNY-00017-1100022



Description :

Nitrile Exam Gloves are reliable, easy-to-wear gloves that are more economical than standard nitrile gloves, and are suited for a wide variety of tasks.

The nitrile gloves deliver enhanced dexterity, thanks to textured fingertips, and are not made with natural rubber latex.

Nitrile Gloves can be an important piece for creating a more efficient, responsible workplace.

The gloves are static dissipative in use. They come in an innovative dispenser box, which incorporates two separate openings.

The first, smaller opening is used when the box is full and helps control dispensing to lessen waste from dropped gloves or multiple dispensing.

When the box is half empty, the second, larger opening allows easier access to the gloves. They have an 1.5 for pinholes and can even be recycled with the Kimberly-Clark Professional RightCycle program, the first large-scale recycling program for nonhazardous lab and cleanroom waste.

Exceptional value: 250% more gloves per case (versus standard nitrile gloves)

Waste less with the one-at-a-time dispenser box; Element Nitrile Gloves are recyclable through the Right Cycle program.

Not made with natural rubber latex, are powder -free and feature beaded cuffs

Superior quality.

Technical Specification :

Professional Nitrile Gloves



Laboratory Instrument India