



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Round Bouffant Caps

**Product Code :**  
LBNY-00017-1100035



**Description :**

Strong, lightweight and breathable  
Latex Free  
Hand-made, flat round style, dual elastic  
Single use  
Comfortable, non-binding elastic band holds the cap securely on the head.  
10-12 gm2 spunbonded polypropylene (SBPP)  
Round bouffants  
Color: white.

**Technical Specification :**

Round Bouffant Caps



**Laboratory Instrument India**