



Email : sales@laboratoryinstrumentindia.com

Product Name :
25 gsm Disposable Isolation Gown

Product Code :
LBNY-00017-1100036



Description :

These disposable polypropylene gowns are used to protect the wearer from contact with potentially infectious liquid and solid material.

Disposable spunbonded polypropylene isolation gown.

Greener Choice Claims

Environmental benefits include:

Non-Toxic.

Free-of a substance that causes environmental harm.

Technical Specification :

25 gsm Disposable Isolation Gown



Laboratory Instrument India