



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Unlined Pigskin Leather Drivers Gloves

**Product Code :**  
LBNY-00017-1100041



**Description :**

Suitable for construction, machine operation, utility work and many general applications  
Gloves return to their original pliability even after being soaked in water.  
Elastic wrist provides a sure fit  
Gunn-cut style with straight thumb  
Unlined.

**Technical Specification :**

Unlined Pigskin Leather Drivers Gloves



**Laboratory Instrument India**