



Email : sales@laboratoryinstrumentindia.com

Product Name :
Low-Lint Nylon Mens Inspection Gloves

Product Code :
LBNY-00017-1100042



Description :

Moisture-absorbent gloves are constructed of 40-denier nylon tricot fabric and can be worn alone, or used as a glove liner.
Easy size identification by color-coded hem on rolled cuff
Can be worn alone or used as a liner
40-denier tricot fabric in a two-piece construction
Moisture absorbent.

Technical Specification :

Low-Lint Nylon Mens Inspection Gloves



Laboratory Instrument India