



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lightweight Elbow-Length Cotton Lisle Inspection Gloves

Product Code :
LBNY-00017-1100044



Description :

Premium quality, cotton lisle inspection gloves
Designed to minimize contamination of products, with an extended cuff for additional protection.
100% cotton fabric
Two-piece, reversible pattern
Effective as liners under disposable and unsupported liquid-proof gloves
Unhemmed.

Technical Specification :

Lightweight Elbow-Length Cotton Lisle Inspection Gloves



Laboratory Instrument India