



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sensor and Transducer Training Kits

Product Code :
LBNY-0001-970001



Description :

Sensor and Transducer Training Kits, Sensor and Transducer Training Kits Equipments, Sensor and Transducer Training Kits Tools, Sensor and Transducer Training Kits Tool Kits, Sensor and Transducer Training Kits Manufacturers, Sensor and Transducer Training Kits Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara

Technical Specification :

Sensor and Transducer Training Kits
To study of humidity sensor transducer
Humidity transducer capacitive for 90% non condensed range 0 to 100%
RH humidity chamber
3.5 digit display voltmeter thru 20 turn pot/preset.
Sockets at different points for observations
Built in IC regulated power supply
Interconnection : 2/4mm banana sockets
220V \pm 10%, 50Hz mains operated
Enclosed in an attractive ABS plastic cabinet with cover

Standard Accessories : User's Manual with patch cords.



Laboratory Instrument India