



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flock-Lined Neoprene/Latex Double-Dipped Gloves

Product Code :
LBNY-00017-1100050



Description :

Protection with comfort
Made of blue neoprene over yellow latex
Raised diamond grip on palm and fingers for good wet or dry grip
Flock lining to absorb perspiration; permits easy on and off
Ideal for general housekeeping, maintenance applications, food processing and food service applications, as well as industrial use.
Straight cuff
Double-dipped construction provides increased protection, versatile chemical resistance properties, and comfort.
Length: 12 in. (30.5 cm)
Thickness: 28 mil.

Technical Specification :

Flock-Lined Neoprene/Latex Double-Dipped Gloves



Laboratory Instrument India

