



Email : sales@laboratoryinstrumentindia.com

Product Name :
Acrylic Glove Box Holder

Product Code :
LBNY-00017-1100051



Description :

Heat-bent walls simplify construction and limit accidental breakage
Triple and quad holders allow glove boxes to be loaded from either the top or side
Clean line storage solutions perfect for laboratory organization, with heat-bent walls to simplify construction and prevent accidental breakage.
Holder offers versatility while also providing space-saving storage
Fits standard glove boxes, accommodating depths up to 3.75 in.
Can be mounted vertically or horizontally, hardware is included.

Technical Specification :

Acrylic Glove Box Holder



Laboratory Instrument India