



Email : sales@laboratoryinstrumentindia.com

Product Name :
Green 15 mil Unlined Chemical Resistant Nitrile Gloves

Product Code :
LBNY-00017-1100055



Description :

Excellent grip with wet and dry applications
Embossed raised-diamond pattern on the palm and fingers
Acceptable for use in both food and pharmaceutical applications.
Color: green.
Size is embossed on each glove.

Technical Specification :

Green 15 mil Unlined Chemical Resistant Nitrile Gloves



Laboratory Instrument India