



Email : sales@laboratoryinstrumentindia.com

Product Name :
Double-Dipped PVC Gloves, Gauntlet Style

Product Code :
LBNY-00017-1100057



Description :

Provides enhanced cut resistance and excellent wet or slippery grip
Enhanced cut resistance
Excellent wet or slippery grip
Interlock lining provides great dexterity
Gauntlet length, covers forearm
Available with sandy or rough acid grip finish
Jersey lining provides comfort and moderate insulation.

Technical Specification :

Double-Dipped PVC Gloves, Gauntlet Style



Laboratory Instrument India