



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Double-Dipped PVC Gloves, Gauntlet Style

**Product Code :**  
LBNY-00017-1100057



**Description :**

Provides enhanced cut resistance and excellent wet or slippery grip  
Enhanced cut resistance  
Excellent wet or slippery grip  
Interlock lining provides great dexterity  
Gauntlet length, covers forearm  
Available with sandy or rough acid grip finish  
Jersey lining provides comfort and moderate insulation.

**Technical Specification :**

Double-Dipped PVC Gloves, Gauntlet Style



**Laboratory Instrument India**