



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Unisex Lab Coats With Knit Cuffs

**Product Code :**  
LBNY-00017-1100059



**Description :**

Knit wrist cuffs provide comfortable and secure fit  
Made of 65/35 Polycotton  
Front snap closure  
Knit cuffs  
One outside left breast pocket, two outside lower pockets  
Standard pointed collar  
White.

**Technical Specification :**

Unisex Lab Coats With Knit Cuffs



**Laboratory Instrument India**