



Email : sales@laboratoryinstrumentindia.com

Product Name :
Disposable SMS Laboratory Jackets

Product Code :
LBNY-00017-1100066



Description :

Premium protective apparel is available in both lab coats and hip length jackets. Made from high quality, ultra-soft SMS fabric designed to keep you protected and comfortable.

The Disposable SMS Lab Jackets with a knitted collar and cuffs. Designed to protect healthcare professionals, lab technicians and other users from fluids and cross-contamination.

Anti-static coating

High quality front snap button closure

Fluid-resistant material meets AAMI L2 standards

42 GSM high quality, breathable SMS material

3-pocket design

Comfortable knit collar and cuffs

Convenient hip length.

Technical Specification :

Disposable SMS Laboratory Jackets



Laboratory Instrument India

