



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Poly/Cotton Bib Apron

**Product Code :**  
LBNY-00017-1100068



**Description :**

Bib apron protects clothing from particles.  
Tubular-braid waist ties  
White.  
Ties in front for easy securing.  
7.5 oz., 65% polyester/35% cotton twill Soil-release finish.

**Technical Specification :**

Poly/Cotton Bib Apron



**Laboratory Instrument India**