



Email : sales@laboratoryinstrumentindia.com

Product Name :
Synthetic Blend String Knit Gloves

Product Code :
LBNY-00017-1100071



Description :

Ideal for use in warehousing, parts handling, assembly, and food processing applications.
Can also be used as glove liners under unprotected work gloves.
Reversible design allows replacement of one glove at a time, rather than pairs
Knit pattern allows breathability with less hand fatigue
Color-coded hem denotes size.

Technical Specification :

Synthetic Blend String Knit Gloves



Laboratory Instrument India