



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Smooth Disposable CPE Gown

**Product Code :**  
LBNY-00017-1100073



**Description :**

Provides light-weight, impervious protection against bio-hazards, with non-linting construction for greater comfort.  
Superior physical properties over imports  
In compliance with FDA requirements and tested in accordance, including breaking load, tensile strength, and elongation  
Successfully tested for liquid and aerosol blood resistance, synthetic blood penetration and resistance, bacteriophage penetration resistance  
Redirects stress to the film material not the thermal seals to minimize breakage potential  
Provides liquid and aerosol blood resistance  
Full-length open back affords better comfort and shoulder coverage  
Lightweight, non-linting construction.  
XL size with easy, tear-away feature for quick gown removal and safe disposal  
Contoured thumb loops secure gown sleeve for added protection.

**Technical Specification :**

Smooth Disposable CPE Gown

---

# Laboratory Instrument India

