



Email : sales@laboratoryinstrumentindia.com

Product Name :
Smooth Disposable CPE Gown

Product Code :
LBNY-00017-1100073



Description :

Provides light-weight, impervious protection against bio-hazards, with non-linting construction for greater comfort.
Superior physical properties over imports
In compliance with FDA requirements and tested in accordance, including breaking load, tensile strength, and elongation
Successfully tested for liquid and aerosol blood resistance, synthetic blood penetration and resistance, bacteriophage penetration resistance
Redirects stress to the film material not the thermal seals to minimize breakage potential
Provides liquid and aerosol blood resistance
Full-length open back affords better comfort and shoulder coverage
Lightweight, non-linting construction.
XL size with easy, tear-away feature for quick gown removal and safe disposal
Contoured thumb loops secure gown sleeve for added protection.

Technical Specification :

Smooth Disposable CPE Gown

Laboratory Instrument India

