



Email : sales@laboratoryinstrumentindia.com

Product Name :
Traffic Signal Control

Product Code :
LBNY-0001-970004



Description :

Traffic Signal Control

Technical Specification :

4 Sequence switch/ 4delay switch
Sequence 1: LED Testing
Sequence 2: One way traffic
Sequence 3 & 4: Two way traffic

Function of LEDs in the trainer.

Green – Start,

Red – Stop

Yellow – Warning To Stop

Blue – Pedestrian

Multi Color– Free Left

Built in Regulated Power Supply: DC +5V/300 mA.

Power Supply:

Voltage range: AC 100V – 230V.

Frequency Range: 50 – 60Hz.

Housing: It is mounted in an elegant ABS Plastic cabinet for better viewing and portability.



Laboratory Instrument India