




Email : sales@laboratoryinstrumentindia.com

Product Name :
Professional Breathable Particle Protection Bouffant Cap

Product Code :
LBNY-00017-1100077



Description :

Made with breathable Microforce barrier fabric
Features serged seams and an elastic opening
Elastic is not made with natural rubber latex (reducing risk for individuals with a Type 1 latex allergy)
Lightweight and breathable 24 in. bouffant caps are best suited for use in cleanroom applications.
Greener Choice Claims
Passes NFPA 99 for antistatic (anti stat) materials
Recyclable via The Right ycle Program by Kimberly-Clark Professional
Environmental benefits include:
Recyclable.

Technical Specification :

Professional Breathable Particle Protection Bouffant Cap



Laboratory Instrument India