




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Professional Breathable Particle Protection Bouffant Cap

**Product Code :**  
LBNY-00017-1100077



**Description :**

Made with breathable Microforce barrier fabric  
Features serged seams and an elastic opening  
Elastic is not made with natural rubber latex (reducing risk for individuals with a Type 1 latex allergy)  
Lightweight and breathable 24 in. bouffant caps are best suited for use in cleanroom applications.  
Greener Choice Claims  
Passes NFPA 99 for antistatic (anti stat) materials  
Recyclable via The Right ycle Program by Kimberly-Clark Professional  
Environmental benefits include:  
Recyclable.

**Technical Specification :**

Professional Breathable Particle Protection Bouffant Cap



**Laboratory Instrument India**