



Email : sales@laboratoryinstrumentindia.com

Product Name :
Professional Sterling 9.5 in. Nitrile Gloves

Product Code :
LBNY-00017-1100078



Description :

Sterling gloves combine the sensitivity of latex with the protection and strength of nitrile.
Typical uses include: Healthcare, Laboratory, Cleaning, Maintenance, Electronic Assembly and Inspection.
Standard 9.5 in. length
Textured fingertips
Beaded cuff
Powder free
Non-sterile
Static dissipative
Ambidextrous.

Technical Specification :

Professional Sterling 9.5 in. Nitrile Gloves



Laboratory Instrument India