



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Professional Sterling 9.5 in. Nitrile Gloves

**Product Code :**  
LBNY-00017-1100078



**Description :**

Sterling gloves combine the sensitivity of latex with the protection and strength of nitrile.  
Typical uses include: Healthcare, Laboratory, Cleaning, Maintenance, Electronic Assembly and Inspection.  
Standard 9.5 in. length  
Textured fingertips  
Beaded cuff  
Powder free  
Non-sterile  
Static dissipative  
Ambidextrous.

**Technical Specification :**

Professional Sterling 9.5 in. Nitrile Gloves



**Laboratory Instrument India**