



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vinyl Closed-Back Aprons, Waist ties

**Product Code :**  
LBNY-00017-1100082



**Description :**

With heat-sealed elastic cuffs  
6 mil thick aprons have electronically welded binding strip around neck edge only; other edges are unhemmed  
Blue  
Neck ties, 15 in. long (38cm); waist ties, 29 in. long (74cm)  
8 mil thick aprons have electronically welded vinyl binding strip that covers all edges for smoother finish.

**Technical Specification :**

Vinyl Closed-Back Aprons, Waist ties



**Laboratory Instrument India**