



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Frocks, Clean-Processed and Sterile

**Product Code :**  
LBNY-00017-1100083



**Description :**

Serged seams have multiple interlocking threads that are sewn around the raw edges of garment material to create a strong, stress-resistant seam

Elastic openings for tighter fit at wrist

Bound neck for lower particle shedding.

Controlled environment processed to protect against particulates and other contaminants. Garments are gamma sterilized to SAL of 10<sup>-6</sup>. Clean-processed garments offer lowest level of particle shedding portfolio.

Front zipper closure for easy donning and doffing

A-line silhouette

Raglan sleeve design for greater range of motion and comfort

Greener Choice Claims

Environmental benefits include:

Recyclable.

**Technical Specification :**

Frocks, Clean-Processed and Sterile



**Laboratory Instrument India**

