



Email : sales@laboratoryinstrumentindia.com

Product Name :

A5 Cleanroom Non-Sterile Hoods

Product Code :

LBNY-00017-1100085



Description :

Non-sterile, white
Triple bagged, vacuum packed
Suitable for Class 5 or higher cleanrooms;
Works with a variety of masks
Clean-processes and individually packaged under cleanroom conditions
Non-sterile cleanroom hoods cover the head, neck, and shoulders to prevent contamination in the environment.
Available with or without ties.
Latex-free
Greener Choice Claims
Environmental benefits include:
Recyclable.

Technical Specification :

A5 Cleanroom Non-Sterile Hoods



Laboratory Instrument India

