



Email : sales@laboratoryinstrumentindia.com

Product Name :
AC Voltage Regulator

Product Code :
LBNY-0002-620020



Description :

AC Voltage Regulator

Technical Specification :

Heavy-duty continuous use transformer.
On/Off switch with indicator lamp, integrated automatic protection circuit.
Analog meter on front of unit measures exact input voltage.
Voltage: 0-400V.
Power: 3000W.
Frequency: 50Hz.

AC Voltage Regulator, AC Voltage Regulator Equipments, AC Voltage Regulator Tools, AC Voltage Regulator Tool Kits, AC Voltage Regulator Manufacturers, AC Voltage Regulator Suppliers from India, China, Kenya



Laboratory Instrument India