



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
A10 Polyethylene Light-Duty Apron

**Product Code :**  
LBNY-00017-1100086



**Description :**

Made of embossed polyethylene for excellent flexibility  
Easy to wipe down and keep clean  
Economical, disposable protection  
Superior resistance to chemicals, oils, fats, and grease.

**Technical Specification :**

A10 Polyethylene Light-Duty Apron



**Laboratory Instrument India**