



Email : sales@laboratoryinstrumentindia.com

Product Name :
Universal Precautions Lab Coats

Product Code :
LBNY-00017-1100087



Description :

Features knit cuffs for increased protection and a traditional collar; Back vent helps keep wearer comfortable and cool

These protective lab coats have 6 snaps in the front and 3 reinforced pockets

Included neck loop means your employees can hang the lab coat wherever it's convenient

The mid-calf length protective coat has knit cuffs, for greater protection, and a traditional collar. Plus, the lab jacket has three reinforced pockets, six snaps and a back vent for increased comfort and ventilation. The neck loop makes it convenient to hang, whether you're using a lab coat in a research laboratory or an industrial environment. Lab Coats with Knit Cuffs + Extra Protection from Kimberly-Clark give your team a blend of comfort and protection in environments where the potential of exposure is low to moderate for chemicals or fluids. In fact, the 3-layer SMS fabric, a nonwoven fabric offers exceptional resistance to blood strike-through at various levels of exposure. Kimberly-Clark Professional is a respected name in protective apparel, enabling you to deliver the reliable protection you're required to provide and the peace of mind your team deserves.

Each unisex lab coat is made of a 3-layer (nonwoven) fabric, for wear when there is potential low to moderate exposure to chemicals or fluids.

Technical Specification :

Universal Precautions Lab Coats

Laboratory Instrument India

