



Email : sales@laboratoryinstrumentindia.com

Product Name :
Coveralls, Sterile

Product Code :
LBNY-00017-1100088



Description :

Coverall has elastic openings for tighter fit at wrist and ankle
Bound neck for lower particle shedding
Front zipper closure for easy donning and doffing
Sterile coveralls used in controlled environments to protect against particulates and other contaminants.
Garments are gamma sterilized of 10-6
Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration
Coverall has elastic openings for tighter fit at wrist and ankle
Bound neck for lower particle shedding
Front zipper closure for easy donning and doffing
Sleeve design for greater range of motion and comfort
Full traceability on all sterilized apparel with certificates of sterility available
Garments are gamma sterilized of 10-6
Bound seams are covered with garment fabric to reinforce the seam and to reduce the potential for particle penetration
Sleeve design for greater range of motion and comfort.

Technical Specification :

Coverall, Sterile



Laboratory Instrument India