



Email : sales@laboratoryinstrumentindia.com

Product Name :
Liquid Barrier Gowns

Product Code :
LBNY-00017-1100089



Description :

The material offers enhanced liquid barrier protection– fabric tested for blood and blood-borne pathogen penetration, chemotherapy drug and hormone permeation and the built-in thumb loops help maintain glove and garment interface.

The gowns are low lint to minimize contamination. Liquid Barrier Gown is bulk-packaged and available in a range of sizes to help ensure comfort.

Liquid Barrier Gowns deliver stringent barrier protection and superior contamination control for your facility. This protective gown has plenty of innovative features, including a seamless front design and an easy to don-and-doff design to reduce contamination risk and prevent hazardous liquids from getting on clothing or skin.

Liquid Barrier Gowns/Case

Seamless front and closed back; Easy to don and doff design

Enhanced liquid barrier protection- fabric tested for blood and blood-borne pathogen penetration, chemotherapy drug and hormone permeation to help protect people and processes from liquid splashes

Thumb loops to maintain glove and garment interface; survey, thumb loops are preferred among compounding pharmacists.

Technical Specification :

Liquid Barrier Gowns

Laboratory Instrument India

