



Email : sales@laboratoryinstrumentindia.com

Product Name :
A5 Sterile Cleanroom Apparel

Product Code :
LBNY-00017-1100090



Description :

The Reflex Coverall Design is 7.5 times more rip-resistant than ANSI minimum standards and offers excellent fit and freedom of movement for the wearer. This personal protective apparel is triple-bagged and vacuum packaged, suitable for Class 5 or higher cleanrooms.

A5 Sterile Cleanroom Coveralls have an average bacterial filtration efficiency (BFE) of 96% and are made of 100% polypropylene breathable SMS fabric, which gives the wearer cloth-like comfort.

It is not made with natural rubber latex, silicone, or the BHT preservative and is designed for single-use (disposable PPE). The easy-to-don, individually packaged coveralls feature a mandarin collar, thumb loops and tunneled elastics.

A5 Sterile Cleanroom Coveralls/Case

Reflex Coverall Design offers excellent fit and freedom of movement and 7.5 times more rip-resistant than ANSI minimums

Triple bagged, vacuum packed, suitable for Class 5 or higher cleanrooms.

This protective apparel isn't made with natural rubber latex, BHT preservative or silicone.

Technical Specification :

A5 Sterile Cleanroom Apparel

Laboratory Instrument India

