



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
White Disposable Visitor Lab Coat

**Product Code :**  
LBNY-00017-1100091



**Description :**

The Disposable Protective Lab Coat is anti-static treated and can be used in light duty applications.  
Designed for light duty applications.  
Offers a barrier against non-hazardous particles and liquid splashes.  
Breathable, lightweight material.  
Anti-static treated.  
Knitted cuffs.  
Zip fastener.

**Technical Specification :**

White Disposable Visitor Lab Coat



**Laboratory Instrument India**