



Email : sales@laboratoryinstrumentindia.com

Product Name :
White Disposable Visitor Lab Coat

Product Code :
LBNY-00017-1100091



Description :

The Disposable Protective Lab Coat is anti-static treated and can be used in light duty applications.
Designed for light duty applications.
Offers a barrier against non-hazardous particles and liquid splashes.
Breathable, lightweight material.
Anti-static treated.
Knitted cuffs.
Zip fastener.

Technical Specification :

White Disposable Visitor Lab Coat



Laboratory Instrument India