



Email : sales@laboratoryinstrumentindia.com

Product Name :
White Disposable Visitor Lab Coat

Product Code :
LBNY-00017-1100091



Description :

The Disposable Protective Lab Coat is anti-static treated and can be used in light duty applications. Designed for light duty applications. Offers a barrier against non-hazardous particles and liquid splashes. Breathable, lightweight material. Anti-static treated. Knitted cuffs. Zip fastener.

Technical Specification :

White Disposable Visitor Lab Coat



Laboratory Instrument India