



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
High Performance Lab Coat, White Frost

**Product Code :**  
LBNY-00017-1100094



**Description :**

Made from high-quality 45 g/m2 SMS fabric  
Detailed finishing, 3-pocket design and front snap closure  
The High Performance Lab Coat, White Frost is premium protective apparel made from top-quality materials and designed to provide reliable barrier protection without sacrificing comfort or style.  
Fluid-resistant fabric Level 2 standards  
Each bag contains 10+2 units, for a total of 12 vs. 10 per package from the competition.

**Technical Specification :**

High Performance Lab Coat, White Frost



**Laboratory Instrument India**