



Email : sales@laboratoryinstrumentindia.com

Product Name :
Jacket Lab Coat

Product Code :
LBNY-00017-1100095



Description :

Material Level 2 requirements for impact penetration and hydrostatic pressure. Latex free
Designed to fit properly with knitted cuffs & collar
High-quality front snap closure buttons
High quality ultra-soft & light weight SMS fabric (Weight 45 g/m² for Jacket & Lab Coat)
Procedure SMS Jacket designed to protect health professionals and patients from cross contamination. Single use and disposable gowns made of polypropylene.
Fluid-resistant material meets AAMI Level 2
Available in a wide range of sizes and in fashionable colors
12 units in every pack.

Technical Specification :

Jacket Lab Coat



Laboratory Instrument India