



Email : sales@laboratoryinstrumentindia.com

Product Name :
Multi-Meter with

Product Code :
LBNY-0001-980001



Description :

Multi-Meter, Multi-Meter Equipments, Multi-Meter Tools, Multi-Meter Tool Kits, Multi-Meter Manufacturers, Multi-Meter Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara

Technical Specification :

As a general measurement tool, widely used in the school, laboratory, and other social fields. can be used on troubleshoot a variety of automotive computer circuit rotary switch.
The electronic meter can measure AC/DC voltage, AC/DC current, resistance, frequency, continuity, diode and Data retention, Manual range, etc.
Multimeter overload protection: Used the PTC protection circuit for resistance, frequency measurement double current fuses. High quality, more safety ant Automatic shutdown.
The digital multimeter with large digital display and color-coded LED's for battery quick check.
The multimeter is a small and compact design, comfortable to handle. It is made of superior material, durable and practical.



Laboratory Instrument India