



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Disposable Long Sleeve Cover

**Product Code :**  
LBNY-00017-1100096



**Description :**

Protective sleeves available in sterile or non-sterile  
Low linting, anti-static fabric  
Extra long length sleeve covers are constructed from antistatic, low-linting material.  
Elasticated arms and cuffs  
Universal size  
Overall length: 500 mm.

**Technical Specification :**

Disposable Long Sleeve Cover



**Laboratory Instrument India**