



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hooded Coveralls

Product Code :
LBNY-00017-1100097



Description :

Coveralls offer an optimum fit with full freedom of movement.

Serged and over-taped seams that provide the performance of Type 3/4/5/6 garments, which have been tested to standards against pressurized liquids, heavy liquid aerosols and airborne solid particles. The coveralls provide an effective barrier against many water-based inorganic chemicals in low concentration and particles, as well as oil repellency.

Providing a suitable protective barrier against many low-concentration, water-based, inorganic chemicals and particles, down to one micron in size,

Elastic face, wrists and ankles as well as glued-in waist elastic

Elastic thumb loops keep sleeves in place

Distinct orange taping of the garment seams providing

Chemical protective clothing, Category III, Type 3-B, 4-B, 5-B and 6-B

Protection against infective agents including resistance to penetration by blood and body fluids using synthetic blood and resistance to penetration by blood-borne pathogens using bacteriophage

Self-adhesive chin flap for tight seal of suit to the mask
color contrast

Fabric and seams offer chemical permeation barrier to low concentration water-based inorganic chemicals

Serged and over-taped seams for protection and strength.

Technical Specification :

Hooded Coveralls



Laboratory Instrument India