



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Mens High-Performance White Lab Coats

**Product Code :**  
LBNY-00017-1100098



**Description :**

With five-button closure  
One-piece, lined, notched lapel collar  
Left breast pocket, barracked pencil stall, two lower pockets  
Side vent openings  
5 oz. performance blend poplin  
80% polyester/20% combed cotton blend  
Elasticated arms and cuffs  
Closure: five buttons  
White.  
Protective sleeves available in sterile or non-sterile  
Low linting, anti-static fabric  
Pre-cure durable press finish  
Care: light soil wash  
Universal size  
Overall length: 500 mm

**Technical Specification :**

Mens High-Performance White Lab Coats

---

# Laboratory Instrument India

