



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mens High-Performance White Lab Coats

Product Code :
LBNY-00017-1100098



Description :

With five-button closure
One-piece, lined, notched lapel collar
Left breast pocket, barracked pencil stall, two lower pockets
Side vent openings
5 oz. performance blend poplin
80% polyester/20% combed cotton blend
Elasticated arms and cuffs
Closure: five buttons
White.
Protective sleeves available in sterile or non-sterile
Low linting, anti-static fabric
Pre-cure durable press finish
Care: light soil wash
Universal size
Overall length: 500 mm

Technical Specification :

Mens High-Performance White Lab Coats

Laboratory Instrument India

