



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Disposable Isolation Gown

**Product Code :**  
LBNY-00017-11000101



**Description :**

These disposable polypropylene gowns are used to protect the wearer from contact with potentially infectious liquid and solid material.

Disposable spunbonded polypropylene isolation gown.

Greener Choice Claims

Environmental benefits include:

Non-Toxic

Free-of a substance that causes environmental harm.

**Technical Specification :**

Disposable Isolation Gown



**Laboratory Instrument India**