



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Boot Covers

**Product Code :**  
LBNY-00017-11000104



**Description :**

Coarsely textured interior to help prevent movement inside  
PVC Gripper soles provide excellent barrier to liquid penetration  
Fabric is antistatic and low linting and provides an excellent barrier to dry particulates, nonhazardous liquids and microorganisms.  
18 in. (45.7cm) high Tyvek uppers with elastic leg and ankle openings  
White  
Serged seams.

**Technical Specification :**

Boot Covers



**Laboratory Instrument India**