



Email : sales@laboratoryinstrumentindia.com

Product Name :
Boot Covers

Product Code :
LBNY-00017-11000104



Description :

Coarsely textured interior to help prevent movement inside
PVC Gripper soles provide excellent barrier to liquid penetration
Fabric is antistatic and low linting and provides an excellent barrier to dry particulates, nonhazardous liquids and microorganisms.
18 in. (45.7cm) high Tyvek uppers with elastic leg and ankle openings
White
Serged seams.

Technical Specification :

Boot Covers



Laboratory Instrument India