



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mandiles Aprons

Product Code :
LBNY-00017-11000106



Description :

Aprons help keep clean and protected when oil, grease, and chemicals are present
Polyurethane construction is strong, tough and offers tear and puncture resistance yet remains flexible in cold temperatures
Helps keeps workers clean and protected from oil, grease and chemicals
Aprons resistant to abrasion, chemicals, oils and fats
Durable
Dielectrically sealed neck and waist straps are easy to keep clean and resist tear out while the neck strap features an adjustable non-corrosive, non-absorbent plastic buckle
8 to 9 mil die cut polyurethane material that outperforms neoprene, nitrile, vinyl and other films.
Uses: Food processing, Life Sciences.

Technical Specification :

Mandiles Aprons



Laboratory Instrument India

